

# **Hot Plate A3**

**Quality – Made in Germany** 



## Big hot plate for laboratory and scientific use

#### Features:

- Temperature control between 27 °C and 45 °C
- Homogeneous temperature distribution
- Planar surface for optimal heat transfer to heating blocks, culture dishes and slides
- Easy to clean glass surface
- Short heating up time
- Extremely flat with bigger work surface
- Space saving due to separate control unit

#### Technical data:

- Dimensions: 420 x 300 x 23 mm (w x d x h)
- Weight: 5.5 kg
- Heating up time: < 15 min.
- Power supply: 100 - 240 V AC, 50/60 Hz, 140 W
- Temperature adjustable between 27 °C and 45 °C
- Safety class IPX1 (Hot Plate)
- Safety class IPX0 (Controller)
- Cable length between Controller and Hot Plate 1.5 m

## Scope of supply:

- Hot plate
- · Controller with mains cable

REF | 15668

## Additional equipment:

Retractor frame

REF 15942



Rev. 3\_01/2017