

ahn myLab SM-01

Single Plate Magnetic Stirrer 1500 rpm

Usage - in laboratories, pharmacies and schools, mixing liquid industry


Maximum speed - from 15 to 1500 rpm with an interval of 50 rpm / Software driven gradual acceleration

Time setting - ranging from 1- 99 minutes / continuous operation feature

Improved mixing - pulse mode where the mixing direction changes every 30 seconds

Maintenance free design - water spill, corrosion and chemically resistant / includes no moving parts / high quality stirring bar

User-friendly interface and design - digital display showing the real-time speed and the remaining time / last run memory function / extremely thin (only 16.0 mm thick)/flat surface enables easy wiping and cleaning



Ordering information

| Description | | Cat. No. | Specifications | |
|---|--|------------|-----------------------------|--------------|
| AHN myLab® SM-01 Single Plate Magnetic Stirrer 1500 rpm | | 7-300-00-0 | Max. stirring quantity | 800 mL |
| | | | Recommended stir bar length | 25 mm |
| Optional accessories | | Cat. No. | Protection class | IP65 |
| Magnetic stirring bar - PTFE coated, 25 x 8 mm | | CRS-MAG | Dimensions (WxDxH) | 20x164x16 mm |
| | | | Weight | 400 g |
| Delivery Package | | | | |
| 1 pc. Magnetic Stirrer 1500 rpm | | | | |
| 1 pc. Stir bar (PTFE coated, 25 x 8 mm) | | | | |
| 1 pc. Instruction manual | | | | |
| 1 pc. Warranty card | | | | |

