

# ahn myLab TR-02

## Tube Rotator 30 rpm

**Usage** - ideal for processing of blood samples


**Capacity** - 16 tubes (4 ESR Tubes and 12 blood collection tubes) with the maximum operational weight capacity of 1.5 kgs

**Fixed speed** - 30 rpm with an accuracy of  $\pm 2$  rpm

**Optimal mixing** - thanks to gentle circular rotations on the fixed plate at the angle of 38 degrees

**Robust and space-saving design** - robust metal construction with rubber feet for stable operation / compact design requires minimum of bench space

**Easy maintenance and adjustment** - the disc plate can be easily removed for cleaning and replacing



## Ordering information

Description	Cat. No.
AHN myLab® TR-02 Tube Rotator 30 rpm	7-210-00-0

Specifications	
Run time	Continuous / infinite
Noise level	Noiseless operation
Max. loading capacity	1.5 kg
Protection class	IP 21
Dimensions (WxDxH)	278x183x252 mm
Weight	1.7 kg

Delivery Package
1 pc. Tube Rotator 30 rpm
1 pc. Power supply adaptor
1 pc. Instruction manual
1 pc. Warranty card

